



REHASENSE



LIGHTWEIGHT
FUNCTIONALITY



SERVER

LIGHTWEIGHT ROLLATOR

The stylish **Server** rollator has a unique triangle profile ensuring a steady frame design. It is among the lightest and easiest to turn rollators on the market in the 150 kg user category.

SERVER

LIGHTWEIGHT ROLLATOR



REHASENSE

FOLDABLE AND EASY TO CARRY

Server is easy to fold and stays folded with an easy-to-use locking system.

Ergonomically designed handle grips make **Server** easy to transport, if needed.

STEADY AND DURABLE

The frame and seat are tested and approved for the maximum user weight of 150 kg. Light and reliable brakes with TPE-covered unbreakable wheels, create a steady and safe driving feeling.

HEIGHT ADJUSTABLE

Server push handles height adjustable to 10 positions, from 63 cm to 102 cm. User comfort is ensured with three different seat heights - 62, 55 and 50 cm - as well as optimal frame design and ergonomic push handles.

ACCESSORIES

Server can be fitted with a broad range of accessories from Rehasense.



Smooth adjustable handle height and lockable handbrakes.



Seat heights 55 and 62 cm.



Maintenance free brakes and TPE-covered unbreakable wheels.



Unique triangle aluminium profile ensures the durability up to 150 kg.



Stays firmly closed with the lock and is easy to open with the push button.



Practical easy-to-remove shopping bag with special designed easy-to-go clips.

TECHNICAL INFORMATION

Colors:



Max. user weight: 150 kg

Item no:

Large
SRWRL60UK (Wine Red)
SRSGL60UK (Sober Grey)

Medium
SRWRM55UK (Wine Red)
SRSGM55UK (Sober Grey)

Measurements:

(width x depth x height)

Large: 60 x 66 x 74–102 cm
(users height 150–200 cm)

Medium: 60 x 66 x 66–86 cm
(users height 135–170 cm)

Seat width: 46 cm

Transportation
measurements folded:
(depth x height x width)

Large: 66 x 79 x 22,5 cm
Medium: 66 x 74 x 22,5 cm

Weight:

Large: 7,0 kg (+ basket 340 g)

Medium: 6,8 kg (+ basket 370 g)

Shopping bag:
(width x depth x height)
32 x 15 x 28 cm
Max. load 5 kg.

Frame material: aluminium



Server is CE-marked and fulfils the requirements of EU directive 93/42 EEC for medical devices.

Copyright © Rehasense® 2019

www.rehasense.com

Rehasense UK
Mobility House
Marsh Lane Industrial Park
Hayle, Cornwall
TR27 5JR

info@rehasense.com